Getting Data

You can easily obtain spectra from your Tracer using the S1PXRF software. This enables you to see your spectrum as it is gathered, and do both qualitative and quantitative analysis.

This guide will take you step by step, button click by button click, through using S1PXRF and your Tracer



Getting Data

Setting up for Analysis





This is the home screen after launching S1PXRF. At this point, you are not connected to the instrument



KTI (.PDZ 40.00kV, 12.50µA)

22.08	25.76	29.43	33.11		36.79) '	4
-15.0C	Amb:101.4F	Ra	w:2,224.68	1	/alid:2,147	7.68	1.
9	Shift:0.00000		TRA	CERTU	ırbo		11.
				Pz †0	P: ↓)	9:16 AM 7/14/2013	





Preview...'

22.08	25.76	29.43	33.11	36.79	4
-15.0C	Amb:101.4F	Raw:2	,224.68	Valid:2,147.68	1.
Sh	nift: 0. 00000		TRA	CERTurbo	1.
			▲ [┺ ĦÎ ╦ ♠) 9:17 AM 7/14/201	3





mana mana			were and the second sec		
ah aa					
22.08	25.76	29.43	33.11	36.79	4
15.0C	Amb:101.4F	Raw:2,3	224.68	Valid:2,147.68	1.
Sh	nift: 0.00000		TRACE	RTurbo	1.
				9:17 A	M

7/14/2013







KTI (.PDZ 40.00kV, 12.50µA)

mand	monorman	and a state of the	www.www.		
2 ² .08	25.76	29.43	33.11	36.79	4
-15.0C	Amb:101.4F	Raw:2,	224.68	Valid:2,147.68	
Sh	nift:0.00000		TRACE	RTurbo	//
			- 🕞	18 A 1 · · · · · · · · · · · · · · · · · · ·	M 013

71

S1PXRF - 14-Jul-2013 09:18:21



X



KTI (.PDZ 40.00kV, 12.50µA)

not be found	×
n the instrument.	
Ignore	

			and the second sec	·	
2 ² .08	25.76	29.43	33.11	36.79	4
-15.0C	Amb:101.4F	Raw:2,	224.68	Valid:2,147.68	
Sh	ift:0.00000		TRACE	RTurbo	1.
			A 🛓	・ 10 記 小 9:18 A	M

_



S1PXRF - 14-Jul-2013 09:19:33





71



m Port 🛛 📉 🗙	
ОК	
Cancel	



m Port	
OK	
Cancel	

71.

S1PXRF - 14-Jul-2013 09:20:37

Back Scatter, when selected, provides a safety feature - your instrument will shut off automatically if X-rays are not all accounted for. However, this can also prevent measurements for small objects. You can temporarily turn it off for some measurements

7/14/2013

PC Trigg	ler		
Do	one		
Back Sc	attei		
nation			

3 09.21.48	_		×
t 14.00			
hotos	Ζ 40.00ΚΥ, ΙΖ.50μΑΙ		
nnle			
55.255			
OK ler			
Actual High Voltage (kV)			
0.05			
Actual Anode Current (μΑ)			
0.05			
	Lovy Filtor (1 mil Al 1 mil Ti)		
I - rei	10W FILLER (1 min AI, 1 min II)		_
	2 - INO FIIter		. 1
3 - Red Fi	Iter (1 mil Al, 1 mil Ti, 1 mil	Cu)
🔍 4 - Green I	Filter (1 mil Al, 1 mil TI, 6 mi	I C	u)
22	5 - Blue Filter (1 mil Ti)		
-15.0C Amb:10 Shift-0.00000	1.4F Raw:2,224.68 Valid:2,147. TBACEBTurbo	68	
		9:21 A	M.
		7/14/2	013

7/14/2013

Getting Data

Step 2 Collecting Data

KTI 40.00kV, 14.80µA (.PDZ 40.00kV, 12.50µA)

Click the '<>' button to expand your window to focus on specific elements. You can also press the Up and Down buttons on your keyboard. You can also click (or on tablets touch) and spread your spectra as well

16.27	16.88	17.48	18.08	
-15.0C	Amb:100.8F	Raw:13,067.02	Valid:12,601.06	1.
Shi	ift:0.00000	TRA	CERTurbo	11.
			P:28 AM 7/14/201	3

KTI 40.00kV, 14.80µA (.PDZ 40.00kV, 12.50µA)

16.27	16.88	17.48	18.08	
:-15.0C	Amb:100.8F	Raw:13,067.02	Valid:12,601.06	
Shi	ft:0.00000	TR/	ACERTurbo	
			1 </th <th>M 13</th>	M 13

16.27	16.88	17.48	18.08	
-15.0C	Amb:100.8F	Raw:13,067.02	Valid:12,601.06	/_
Shift:0.00000		TRACERTurbo		/_



S1PXRF - 14-Jul-2013 09:30:41





Х

S1PXRF - 14-Jul-2013 09:30:59

71.



7/14/2013



Always select to save as a 'PDZ' file - this is the file format for other Bruker software.

'TXT' and 'CSV' can be used in Excel and other programs - but this is the raw data, it may only be useful to advanced users

Make sure Autosave is selected

After you have selected settings, press the 'OK' button

r Tm Yb Lu

33.22 22.15 29.53 36.92 25.84 Det:-15.0C Amb:100.8F Raw:13,067.02 Valid:12,601.06 Shift:0.00000 TRACERTurbo 7/14/2013











33.22 22.15 36.91 25.84 29.53

Amb:102.1F Raw:13,072.34 Valid:12,585.11 Shift:0.00000 TRACERTurbo 9:34 AM **†1) 🏝**

7/14/2013







71.

S1PXRF - 14-Jul-2013 09:34:26







				· · · · · · · · · · · · · · · · · · ·	
22.15	2.15 25.84		33.22	36.91	
15.0C	Amb:102.1F	Raw	:13,072.34	Valid:12,585.11	1.
Sh	nift:0.00000		TRA	CERTurbo	1.
				9:34	AM

7/14/2013

X







7/14/2013

Getting Data

Step 3 Navigating and Analyzing Spectra







X

















\sim	<u> </u>	\sim		-		- man				-		~~~
1	6.03	}	7	.23			8.44				9.64	
-15.0C		Amb:1	02.8F		Raw:9	,731.71		٧	alid:	9,439	9.02	1.
	Shift:	0.00000					TRAC	ERTu	rbo			11.
							-	3 †	*	())	9:39 AM	





The next element at first glance seems to be Titanium - but the peaks do not match

\sim	\sim			and the second s			
6	5.03	7	.23	8.4	4 '	9.	64
-15.0C	Amb:1	02.8F	Raw:	9,731.71	Vali	d:9,439.02	1.
Sł	nift:0.00000			TBA	ACERTurbo)	11.
					Pa †0 *	9:3 <u>9</u> :3	9 AM 4/2013















S1PXRF - 14-Jul-2013 09:42:12





X











Getting Data

Step 4 Quantifying Data







22.14	25.84	29.53	33.22	36.91	4
-15.5C	Amb:94.7F	Raw:41,	871.66	Valid: 38,009.63	1.
Shil	ft:0.00000		TRACE	RTurbo	1.
			-	; †⊡ ‱ ()	AM /2013

 \mathbf{X}



22.14	25.84	29.53	33.22	36.91	- 4
-15.5C	Amb:94.7F	Raw:41	,871.66	Valid:38,009.63	/_
	Shift:0.00000		TRACE	RTurbo	

7/14/2013

 \mathbf{X}


71. Setup ID DownLoad Timed Tube Help File real a 🕑 <u>S</u>tar 21 **Read File** PeakID T3S2430 Emperical Cals T3S2430 ▶ ~ C Search T3S2430 € <u>B</u> <u>M</u> (→) << New folder Organize 🔻 Date modified Name 쑦 Favorites 📃 Desktop AL1.CFZ 1/11/2013 11:44 AM 👃 Downloads AL1.SRZ 1/11/2013 11:43 AM 🙊 Photo Stream ALV1.SRZ 1/11/2013 11:44 AM 🔚 Recent places ALV11.CFZ 1/11/2013 11:44 AM ALV12.CFZ 1/11/2013 11:44 AM 📄 Libraries AU1.CFZ 1/11/2013 11:43 AM Documents AU1.SRZ 1/11/2013 11:43 AN 🎝 Music CO1.CFZ 1/11/2013 11:43 AM Pictures CO1.SRZ 1/11/2013 11:43 AM 😽 Videos 🔼 CU1.CFZ 1/11/2013 11:43 AM CU1.SRZ 1/11/2013 11:43 AM Normal Computer CU2.CFZ 1/11/2013 11:43 AM BOOTCAMP (C:) FE1.CFZ 1/11/2013 11:43 AN 👝 Lincoln (D:) FE1.SRZ FE2.CFZ 1/11/2013 📬 Network NI1.CFZ <u>A</u>= 1/11/201 < Coef File name: FE1.CFZ ¥ Open LSec:3.74 Chan:1903.929 Kev:38.2417 Cnt:0 Det: Vac:715 evch:20.0857

71

Po

S1PXRF - 14-Jul-2013 09:47:05

		X	DZ 40.00kV, 14.80µA)
0 E	Emperical Cals	<u>م</u>	Select Elements
	== ▼		Al Si P S Cl Ar Mn Fe Co Ni Cu Zn Ga Ge As Se Br Kr Tc Ru Rh Pd Ag Cd In Sn Sb Te I Xe Be Os Ir Pt Au Ha TL Ph Bi Po At Bn
Л Л	CFZ File		Nd Pm Sm Eu Gd Tb Dy Ho Er Tm Yb Lu
1	SRZ File		U Np Pu Am Cm
1	CFZ File		
1	CFZ File		
1	CFZ File		
1	SRZ File		
1	CFZ File		
1	SRZ File		
1	CFZ File		
1	SRZ File		
1	CFZ File		
1	CFZ File		

Since Iron is the dominant element in this spectra, we will use the FE1.CFZ file

-	Cancel					T	
		29	.53	33.22	36.91	· · ·	4
-15.5C	Amb:S	94.7F	Raw: 41,871.0	66	Valid:38,009	.63	h.
	Shift:0.00000			TRACER	Turbo		k
				▲ ►	†© * 1 •>	9:47 AM 7/14/2013	

 $-\Box$ \times



22.14	25.84	29.53	33.22	36.91	∣ ′ 4
-15.5C	Amb:94.7F	Raw:41,	,871.66	Valid: 38,00	9.63 //
Shi	ift:0.00000		TRACI	ERTurbo	1.
			- F	; †î *i •)	9:47 AM 7/14/2013

 \mathbf{X}



5	Result Table	×						
Copy Save	-Inte							
Hotel Glass Y	ellow Filter 07-14-2013 09:47:25	Concentration (FE						
TiKa1	0.0020568686	0.086126619						
V Ka1	0.0026004382	0.032601909						
CrKa1	0.2237540267	22.725950990						
MnKa1	0.0465116729	1.749850070						
FeKa1	0.5244142514	64.936527879						
CoKa1	0.0376943989	0.531770488						
NiKa1	0.0356112622	2.112258123						
CuKa1	0.0028178025	0.236141518						
W Lb1	0.0003884333	0.120002285						
MoKa1	0.1241508452	3.694419234						
	In the first co includ	lumn, we led in the	have the elements calibration					
NiKa1 CoKa1 U Ka1 U Ka1 TiKa1 CuKa1W Lb1								
3.69	7.38	11.07	14.76 18.45					
it:0	LSec:3.74	Chan:1903.929	Kev:38.2417 Det:					
evch:2	.0.0857		Vac:715					
		L BASSIN						

KTI 40.00kV, 14.80µA (.PDZ 40.00kV, 14.80µA)



7/14/2013

_

22.14	25.84	29.53	33.22	36.91	- 4
Det:-15.5C	Amb:94.7F	Raw:41	,871.66	Valid: 38,009.63	1.
Shi	ift:0.00000		TRACE	RTurbo	1.
			· 🛌	9:47 A	M



				Select E	lements	
Conv. Save lat	Result Table	×	H Li Be Na Mo	TiKa1 · 22 (K)	LOM <u>He</u> Z+ <u>BCNDFNe</u> AlSiPSCIAr	
Copy Save-Int				Y Zr Nh Ma Ta Bu Bh	<u>Ni Lu Zn Galge As Se Br Kr</u> Pd Ag Cd In Sn Sh Tell Xe	
Hotel Glass Yellow	Pilter 07-14-2013 09:47:25	Concentration (FE	Cs Ba	La Hf Ta W Be Os Ir	Pt Au Ha TI Pb Bi Po At Bn	
	0.0020568686	0.086126619	Fr Ba			
	0.0026004382	0.032601909				
	0.2237340267	1 749950070			Amjumj	
	0.0460116723	1.743630070 C4.020E27070				
	0.0244142014	0 521 770 400				
	0.0376343363	2 112250122				
	0.0306112622	0.2261/1510				
	0.0020170020	0.230141310				
W LDI	0.1013004333	0.120002283				
	The	second column ir used to estimate	elemental weig	atical param ght percents	eters	
	NiKa1 ÇeKa1	Mo]	(a1			
Mu V Kai TiKai	Kal CuKal¥ Lb1			· • · · · · · · · · · · · · · · · · · ·		
3.69	7.38	11.07 14.76	18.45 22.14	25.84	29.53 33.22	36.91
nt:0	LSec:3.74	Chan:1903.929 Kev:38.	2417 Det:-15.5C	Amb:94.7F	Raw:41,871.66	Valid: 38,009.63
evch:20.0	857	Vac:715		Shift:0.00000	TRACER	l urbo
(▲ <mark> </mark>	#1 ** (*) 9:4 7/14

KTI 40.00kV, 14.80µA (.PDZ 40.00kV, 14.80µA)

_ 🗆 🗡



5	Result Table	×	
Copy Save-Inte			
Hotel Glass Yellow Filte	07-14-2013 09:47:25	Concentration (FE*	
TiKa1	0.0020568686	0.086126619	
V Ka1	0.0026004382	0.032601909	
CrKa1	0.2237540267	22.725950990	
MnKa1	0.0465116729	1.749850070	
FeKa1	0.5244142514	64.936527879	
CoKa1	0.0376943989	0.531770488	
NiKa1	0.0356112622	2.112258123	
CuKa1	0.0028178025	0.236141518	
W Lb1	0.0003884333	0.120002285	
MoKa1	0.1241508452	3.694419234	
U Kal	NiKa1 Na1	The thir cas Molybd othe	d column includes se, we have ~23% (enum, ~2% Nickel er smaller concentr
TiKa1	CuKa1V Lb1		\dots
3.69	7.38	11.07	14.76 18.45
Cnt:0 LS	ec:3.74	Chan:1903.929	Kev:38.2417 Det:
evch:20.0857	N		

KTI 40.00kV, 14.80µA (.PDZ 40.00kV, 14.80µA)



9:47 AM

7/14/2013

†1) 🏝







KTI 40.00kV, 14.80µA (.PDZ 40.00kV, 14.80µA)

						Se	elec	t E	ler	ner	nts						
H	Be		TiKa Cle	a1 - 2 ar Al	22 [() Ac	K Id	O Z	L	O Z+	M	<u>B</u>	<u>C</u>	N	0	F	<u>He</u> Ne
Na K Rb	<u>Mg</u> Ca Sr	<u>Sc</u> Y	Ti Zr	V Nb	Cr Mo	<u>Mn</u> Tc	Fe Ru	<u>Co</u> Rh	Ni Pd	<u>Cu</u> Ag	J Zn Cd	<u>Al</u> Ga In	<u>Si</u> Ge Sn	P As Sb	<u>S</u> Se Te	<u>CI</u> Br I	<u>Ar</u> Kr Xe
<u>Cs</u> Fr	<u>Ba</u> Ra	<u>La</u> Ac	Hf	Ta Ce	W Pr	Re Nd	Os Pm	lr Sm	Pt Eu	Au Gd	<u>Hg</u> ТЬ	TI Dy	РЬ Ho	Bi Er	Po Tm	Аt Yb	<u>Rn</u> Lu
	<u>0</u> K			Th	Pa	U	Np	Pu	Am	Cm							

7/14/2013

We can select the cells in this window and copy them into Excel, Word, or any other software program

22.14	25.84	29.53	33.22	2 36.9	1 4
-15.5C	Amb:94.7F	Rav	w: 41,871.66	Valid:38,00)9.63 //.
Sh	ift: 0.00000		TR	ACERTurbo	h.
					9:47 AM